

GD4210



Article Code: **GD4210** For 1 Monitor up to 19" for Portegra2 Suspension Arms
GD4210-L For 1 Monitor up to 21/24" for Portegra2 Suspension Arms


CE
 Made in Germany

Portegra2 Mono Monitor Holder - GD4210 Models

Monitor Support for One Flat Screen Display.

The GD4210 can easily hold one monitor up to 21/24". The large curved handlebar protects the monitor and provides a means for quick & easy positioning. Combined with the mobility and safety of the Portegra2 suspension arm, the GD4210 offers the user a worry-free and flexible application.

Type:	Monitor holder with handle and VESA 75/100 adapter for one flat screen display, external cabling	
Weight Load:	GD4210	9.5 kg / 20.9 lbs
(Max.)	GD4210-L	15.4 kg / 33.9 lbs
Weight:	GD4210	2.5 kg / 5.5 lbs
	GD4210-L	2.6 kg / 5.7 lbs
Colour:	RAL 9002	
Test mark:	CE, UL	



* (Dimensions for L version)

MAVIG GmbH
 X-Ray Protection
 & System Solutions

Stahlgruberring 5
 81829 Munich
 Germany

PO Box 82 03 62
 81803 Munich
 Germany

Phone: +49 (0) 89 / 420 96 - 0
 Fax: +49 (0) 89 / 420 96 - 200
 e-Mail: info@mavig.com

www.mavig.com