

PORTEGRA²
Mounting Plate for stationary columns



TS1520 Mounting plate for stationary ceiling columns


CE, IEC 60601-1
 Made in Germany

MAVIG Mounting Plate for Fixpoint Installations

MAVIG has set a standard for practices and hospitals in ceiling support and radiation protective systems.

A solid construction for maximum safety in rooms with hanging substructures. MAVIG's Portegra2 mounting plate is a great solution for an easy installation and reliable use.

- Suitable for all MAVIG Portegra2 stationary ceiling columns
- Solid construction, designed for MAVIG's patented ceiling anchor safety system
- Cable cut-out for extra power, imaging or injector cables
- Recommended for substructures, e.g. Unistrut®
- Dimensions: 40 x 70 cm (16 x 28 in)
- Includes necessary hardware for mounting a Portegra2 ceiling column



MAVIG GmbH
 X-Ray Protection
 & System Solutions

Stahlgruberring 5
 81829 Munich
 Germany

PO Box 82 03 62
 81803 Munich
 Germany

Phone: +49 (0) 89 420 96 - 0
 Fax: +49 (0) 89 420 96 - 200
 e-Mail: info@mavig.com

www.mavig.com